

Options

- **Adjustable Extension Belt** – Available in four different adjustable lengths to extend lifting straps:
 - 12" – 16" (305 mm – 406 mm)
 - 16" – 24" (305 mm – 610 mm)
 - 24" – 40" (610 mm – 1016 mm)
 - 40" – 55" (1016 mm – 1397 mm)
- **Carriages** – Three recommended carriages to mount Multirall onto overhead rail:
 - **D45** – Double hook
 - **S50** – Single hook
 - **S65** – Single hook
- **Easy Switch** – Fits hand control with membrane key for easy and comfortable motor operation
- **Hang-Up Hand Control Hanger** – Provides easy and convenient storage for Multirall hand control
- **Holder for Accessories** – Attaches to *Trolley* for easy storage of *Extension Arm* or *Helping Hand*
- **LikoScale 350** – Scale can weigh patients up to 772 lbs (350 kg)
- **Liko Freespan** – Liko Multirall is compatible with two of the Liko Freespan rail systems, including the:
 - Straight Rail Lift System <product page link>
 - Traverse Lift System <product page link>
- **Liko Support Springs** – Three recommended elastic springs for bounce and gait training:
 - **Short**
 - Weight capacity: 154 lbs (70 kg)
 - Unit length: 8.7" (220 mm)
 - **Short**
 - Weight capacity: 220 lbs (100 kg)
 - Unit length: 8.7" (220 mm)
 - **Long**
 - Weight capacity: 154 lbs (70 kg)
 - Unit length: 11" (280 mm)
- **Parking Panel** – Mounts to wall and provides easy and convenient storage for slingbar, sling and other Multirall accessories
- **Quick-Release Hook** – Attaches to Liko sling bars not already equipped with Quick-Release Hook for safe and easy changing of accessories
- **Reacher arms** – Two recommended arms to aid the connection and disconnection of the lift strap:
 - **Extension Arm** – Compatible with *Carriage S65*
 - **Helping Hand** – Compatible with *Carriage S50*

- **Room-to-Room Strap** – Aids in room-to-room patient transfers with no doorway modifications. Compatible with *Strap Gripper* and available in three different lengths:
 - 46" (1168 mm)
 - 70" (1778 mm)
 - 93" (2362 mm)
- **SlingBar Cover Paddy 30** – Covers sling bar to create padding for patients who risk hitting their heads on sling bar due to involuntary movement. Compatible with:
 - *Universal SlingBar 350*
 - *Universal SlingBar 450*
 - *Universal SlingBar 600*
- **Slings** – A variety of Liko's slings are compatible with Multirall and its sling bars to achieve safe, versatile and comfortable patient suspension and transfer. With over 20 different sling models, sized from XS to bariatric, and many different materials, Liko slings can be used for:
 - General lifting situations
 - Hygiene lifting in connection to toiling or bathing
 - Repositioning
 - Gait training and rehabilitation
 - Bariatric lifting
 - Specialized lifting situations
- **Sling bars** – Seven recommended sling bars for versatile and comfortable patient suspension:
 - **Liko SlingBar Mini 220** – Weight capacity: 450 lbs (205 kg)
 - **Liko Sling Cross-bar 450 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Liko Sling Cross-bar 670 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 350 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 450 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 600** – Weight capacity 660 lbs (300 kg)
 - **Universal TwinBar 670 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
- **Strap Gripper** – Connects to *Room-to-Room Strap* to aid in room transfers
- **Stretchers** – Four recommended stretchers to enable easy horizontal lifting:
 - **FlexoStretch**
 - Weight capacity: 660 lbs (300 kg)

- Unit width (adjustable): 27.6" – 34.4" (700 mm – 950 mm)
 - Unit length:
 - 65" (1650 mm)
 - 73" (1850 mm)
- **LikoStretch Mod 600 IC Lift Stretcher**
 - Weight capacity: 550 lbs (250 kg)
 - Unit width:
 - 28" (711 mm)
 - 32" (813 mm)
- **LikoStretch 1900 Lift Stretcher**
 - Weight capacity: 550 lbs (250 kg)
 - Unit dimensions: 74.80" x 35.79" x 23.35" (1900 mm x 909 mm x 593 mm)
- **OctoStretch Lift Stretcher**
 - Weight capacity: 440 lbs (200 kg)
 - Unit dimensions:
 - 63" x 20" (1600 mm x 500 mm)
 - 67" x 35.4" (1700 mm x 900 mm)
- **Trolley** - Allows for easy storage and transportation of Multirall and its corresponding accessories