

Technical specifications for Liko Sabina 200:

- Weight capacity: 350 lbs (160 kg)
- Lifting distance: 31.3" (795mm)
- Lifting speed (unloaded): 0.9"/s (0.02 m/s)
- Unit dimensions:
 - Minimum: 45.6" x 26.7" x 41.3" (1160 mm x 680 mm x 1050 mm)
 - Maximum: 45.6" x 44" x 68.9" (1160 mm x 1120 mm x 1750 mm)
- Unit weight: 93 lbs (42 kg)
- Turning diameter: 49.2" (1250 mm)
- Battery: 2 x 12V 2.9 Ah

Technical specifications for Liko Sabina II:

- Weight capacity: 440 lbs (200 kg)
- Lifting speed (unloaded):
 - 1.2"/s (0.03 m/s)
 - 1.9"/s (0.05 m/s)
- Unit weight: 100 lbs (45.2 kg)
- Turning diameter: 46.5" (1180 mm)
- Battery: 2 x 12V 2.9 Ah

Options for both Liko Sabina 200 and Liko Sabina II models:

- **Heel Support** – Attaches on footrest for extra support to patient's feet during raising and lowering
- **SeatStrap** – Attaches to *SeatStrap Slingbar* to facilitate first part of raising motion for patient
- **SeatStrap Slingbar** – Compatible with *SeatStrap* for versatile and comfortable patient suspension
- **SideSupport** – Attaches on lift to provide extra lateral stability for patient

Option for Liko Sabina 200 only:

- **Calf strap** – Secures leg position for extra patient stability

Options for Liko Sabina II only:

- **Battery charger, wall-mounted** – Attaches on wall for convenient battery charger storage
- **ComfortSlingBar** – Attaches to lift for versatile and comfortable patient suspension:
 - Lifting distance: 30.9" (785 mm)
 - Liko Sabina II dimensions:
 - Minimum: 41.7" x 27.2" x 41.3" (1060 mm x 690 mm x 1050 mm)
 - Maximum: 41.7" x 44.9" x 68.9" (1060 mm x 1115 mm x 1750 mm)
- **Extra battery** – Lift battery for spare usage

- **SlingBar 350** – Attaches to lift for versatile and comfortable patient suspension:
 - Lifting distance: 32.5" (825 mm)
 - Liko Sabina II dimensions:
 - Minimum: 41.7" x 27.2" x 41.3" (1060 mm x 690 mm x 1050 mm)
 - Maximum: 41.7" x 44.9" x 69.6" (1060 mm x 1115 mm x 1770 mm)
- **Slings** – Three recommended slings are compatible with the Liko Sabina II to achieve safe, versatile and comfortable patient suspension and transfer. Liko offers:
 - ComfortVest sling
 - SafetyVest sling
 - SupportVest sling
- **Table charger housing** – Provides convenient, table-top battery charger storage