

Options

- **Battery Charger** - Lift battery charger for spare usage
- **LikoScale 350** – Scale can weigh patients up to 772 lbs (350 kg)
- **Liko Support Springs** – Three recommended elastic springs for bounce and gait training:
 - **Short**
 - Weight capacity: 154 lbs (70 kg)
 - Unit length: 8.7" (220 mm)
 - **Short**
 - Weight capacity: 220 lbs (100 kg)
 - Unit length: 8.7" (220 mm)
 - **Long**
 - Weight capacity: 154 lbs (70 kg)
 - Unit length: 11" (280 mm)
- **Parking Panel** – Mounts to wall and provides easy and convenient storage for slingbar, sling and other Likorall accessories
- **Quick-Release Hook** – Attaches to Liko sling bars not already equipped with Quick-Release Hook for safe and easy changing of accessories
- **R2R sling bars** – Three recommended sling bars equipped with Liko's R2R technology, which allows room-to-room transfers without doorway modifications:
 - **SlingBar Mini 220 R2R** – Weight capacity: 440 lbs (205 kg)
 - **Universal SlingBar 350 R2R** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 450 R2R** – Weight capacity: 660 lbs (300 kg)
- **SlingBar Cover Paddy 30** – Covers sling bar to create padding for patients who risk hitting their heads on sling bar due to involuntary movement.
Compatible with:
 - *Universal SlingBar 350*
 - *Universal SlingBar 450*
 - *Universal SlingBar 600*
- **Sling bars** – Eight recommended sling bars for versatile and comfortable patient suspension:
 - **Liko SlingBar Mini 220** – Weight capacity: 450 lbs (205 kg)
 - **Liko Sling Cross-bar 450 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Liko Sling Cross-bar 670 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Universal SideBars 450 including bag** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 350 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 450 with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)

- **Universal SlingBar 600 with Quick-Release hook** – Weight capacity 660 lbs (300 kg)
- **Universal SlingBar 670 Twin with Quick-Release hook** – Weight capacity: 660 lbs (300 kg)
- **Stretchers** – Four recommended stretchers to enable easy horizontal lifting:
 - **FlexoStretch**
 - Weight capacity: 660 lbs (300 kg)
 - Unit width (adjustable): 27.6" – 34.4" (700 mm – 950 mm)
 - Unit length:
 - 65" (1650 mm)
 - 73" (1850 mm)
 - **LikoStretch Mod 600 IC Lift Stretcher**
 - Weight capacity: 550 lbs (250 kg)
 - Unit width:
 - 28" (711 mm)
 - 32" (813 mm)
 - **LikoStretch 1900 Lift Stretcher**
 - Weight capacity: 550 lbs (250 kg)
 - Unit dimensions: 74.80" x 35.79" x 23.35" (1900 mm x 909 mm x 593 mm)
 - **OctoStretch Lift Stretcher**
 - Weight capacity: 440 lbs (200 kg)
 - Unit dimensions:
 - 63" x 20" (1600 mm x 500 mm)
 - 67" x 35.4" (1700 mm x 900 mm)